Complex Sameness: Tests to Determine Properties for PLGA Excipients in Long-Acting Formulations

John Garner jg@akinainc.com



3495 Kent Avenue West Lafayette, IN 47906 U.S.A.

Akinainc.com

Akina Overview



- > Research and Development company based in Purdue Research Park
- > Founded by Professor Kinam Park in 2001
- > Two major divisions: PolySciTech & Akinalytics





- Over 650 products available
- Over 375 publications using products



Akinalytics

Advanced analytical polymer assays

- Full characterization of PLGA
- Contracted analysis and research



- > Abbreviated New Drug Application (ANDA) for proposed generic product based on a reference listed drug (RLD)
- > Require demonstration of sameness by *in vitro* methods
 - **Q1: Qualitative sameness.**
 - Q2: Quantitative sameness:
 - **Q3µS:** Microstructural sameness:

Identification of the excipient type Identification of the excipient quantity Identification of the microstructure

PLGA



- Poly(lactide-co-glycolide) (PLGA)
- > Has been used since 1970s: OLD!
- > Properties:
 - Lactide:Glycolide (L:G) Ratio Molecular weight / Distribution End-cap Randomness versus blockiness in glycolide distribution
 - Molecular shape
 - Chirality

Assays must be performed on the final product formulation approved by the FDA



Lactic acid Glycolic acid

Clinical Products using PLGA



Product Name	API*	Туре	Duration	Dose	Approved
Lupron Depot [®]	Leuprolide acetate	Microparticle	1,3,4,6 months	7.5 mg/month	1989, 1996, 1997,2011
Zoladex [®] Depot	Goserelin acetate	Solid implant	1,3 months	3.6 mg/month	1989
Sandostatin [®] LAR	Octreotide acetate	Microparticle	1 month	10-30 mg/month	1998
Atridox®	Doxycycline hyclate	In situ gel*	1 week	42.5 mg/week	1998
Nutropin Depot®	Somatotropin	Microparticle	1 month	13.5 mg/month	1999
Trelstar®	Triptorelin pamoate	Microparticle	1,3,6 months	3.75 mg/month	2000, 2001, 2010
Somatuline [®] Depot	Lanreotide	Microparticle	1 month	60 mg/month	2000
Arestin®	Minocycline HCl	Microparticle	2 weeks	0.5 mg/week	2001
Eligard®	Leuprolide	In situ gel**	1,3,4,6 months	7.5 mg/month	2002
Risperidal [®] Consta [®]	Risperidone	Microparticle	2 weeks	12.5 mg/week	2003
Vivitrol®	Naltrexone	Microparticle	1 month	380 mg/month	2006
Ozurdex®	Dexamethasone	Solid implant	3 months	0.23 mg/month	2009
Propel®	Mometasone furoate	Solid implant	1 month	0.37 mg/month	2011
Bydureon®	Exenatide	Microparticles	1 week	2.0 mg/week	2012
Lupaneta Pack [™]	Leuprolide acetate	Microparticles	3 months	3.75 mg/month	2012
Signifor [®] LAR	Pasireotide	Microparticles	1 month	20~60 mg/month	2014
Zilretta®	Triamcinolone acetoamide	Microparticles	3 month	32 mg/3 months	2017
Sublocade™	Buprenorphine	In situ gel	1 month	100, 300 mg/month	2017
Perseris TM	Risperidone	In situ gel	1 month	90, 120 mg/month	2018

To date, no generic LAI formulation has been approved in USA.





Dissolves API Precipitates polymer

PLGA: Advanced Extraction





Skidmore, Sarah, Justin Hadar, John Garner, Haesun Park, Kinam Park, Yan Wang, and Xiaohui Jiang.

Complex sameness: Separation of mixed poly (lactide-co-glycolide) s based on the lactide: glycolide ratio. Journal of Controlled Release 30 (2019) 174-184.

PLGA: L:G ratio by ¹H NMR







Lactic acid Glycolic acid

 $M_{\rm L} = \frac{P_{\rm L}}{(P_{\rm L} + (P_{\rm G}/2))}$ $M_{\rm G} = \frac{(P_{\rm G}/2)}{(P_{\rm L} + (P_{\rm G}/2))}$

- L:G ratio is critical to the degradation kinetics.
- Higher L:G ratio degrades more slowly

PLGA: Endcap by ¹³C NMR





John Garner, Skidmore, Sarah, Haesun Park, Kinam Park, Stephanie Choi, Yan Wang, A protocol for assay of poly(lactide-co-glycolide) in clinical products. Int. J. Pharm. 495 (2015) 87-92.

PLGA: Blockiness





- L:G react at different rates
- Not uniformly 'random' but reaction method affects distribution
- 'Sameness' or just quality?
- Rcms = GL/GG (Hausberger)
- Rc = GG/GL (Skidmore)
- Hausberger, Angela G., and Patrick P. DeLuca. Characterization of biodegradable poly (D, L-lactide-coglycolide) polymers and microspheres. J. Pharm. Biom. Anal. 13(6) (1995): 747-760.
- Skidmore, Sarah, Justin Hadar, John Garner, Haesun Park, Kinam Park, Yan Wang, and Xiaohui Jiang. Complex sameness: Separation of mixed poly (lactide-co-glycolide) s based on the lactide: glycolide ratio. Journal of Controlled Release 30 (2019) 174-184.

PLGA: Chirality





Scientific Poster Presentation has more details.

- Running HNMR with decoupling (1) separates isotactic lactide (5.17 ppm) from other lactide forms (atactic, syndiotactic, ~5.2-5.3 ppm)
- For PLGA 75:25 or greater, can distinguish P(DL)LG from PLLG. A lower lactide content has a substantial glycolide interference.
 - Zell, Mark T., Brian E. Padden, Amanda J. Paterick, Khalid AM Thakur, Robert T. Kean, Marc A. Hillmyer, and Eric J. Munson. Unambiguous determination of the 13C and 1H NMR stereosequence assignments of polylactide using highresolution solution NMR spectroscopy. Macromolecules 35 (2002) 7700-7707.

PLGA: MW by GPC-ES (External Standard)





34.00 32.00 **Decreasing Mw** 30.00 28.00 26.00 24.00 22.00 20.00 18.00 16.00 14.00 12.00 10.00 8.00 6.00 4.00 2.00 0.00 10.00 11.00 12.00 13.00 17.00 18.00 19.00 20.00 21.00 22.00 23.00 24.00 26.00 14.00 15.00 16.00 25.00

- Size-exclusion chromatography/Gelpermeation chromatography separate based on the hydrodynamic volume
- Solvation between polymer/solvent, interaction with column also controls the retention time.
- Non-representative standard (typically polystyrene) and the lack of standardized methods mean lab-to-lab differences in GPC measured Mw.

Polystyrene Standards (1ml/min THF)

PLGA: MW by GPC-4D





- GPC separation followed by MALLS (Rg, Mw), Viscometry (intrinsic viscosity), inline dynamic light scattering (R_h), and refractive index (concentration).
- Universal calibration, no external standards, no solvation-artifacts
- Light-scattering is better for low-RI solvents (acetone) than high-RI solvents (THF)
- In-depth information about PLGA for determination of molecular shape & branching.

PLGA: Branching Units/Molecule of Sandostatin®

 3.25 ± 0.18

 3.18 ± 0.20

Glu-PLGA (Sandostatin 10)

Glu-PLGA (Sandostatin 28)

 $(0.0980 \,\mathrm{mL/g})$

 $(0.0980 \, \text{mL/g})$





- Sandostatin lists Glucose-PLGA-Glucose without any specific information on branching.
- GPC-4D branched-PLGA method was developed/validated using a series of standards.

 3.07 ± 0.31

 2.75 ± 0.37

• Tested Sandostatin extract: branching average typically ranges between ~ 2.5-4 branching units/molecule.

 2.85 ± 0.41

 2.55 ± 0.49

 $36,676 \pm 1020$

 $39,063 \pm 1561$

 $43,012 \pm 856$

 $46,473 \pm 1248$

 32.19 ± 0.75

 33.97 ± 1.72

 4.77 ± 0.56

 5.11 ± 0.82

 5.67 ± 0.09

 5.83 ± 0.15

Hadar, Justin, Sarah Skidmore, John Garner, Haesun Park, Kinam Park, Yan Wang, Bin Qin, and Xiaohui Jiang. <u>"Characterization of branched poly (lactide-co-glycolide) polymers used in injectable, long-acting formulations.</u>" Journal of Controlled Release (2019).

PLGA: Different PLGAs in Trelstar[®] 22 mg

Semi-solvents for Lacto-selective Solubility

• Simple extraction – NMR only gives average LA:GA of components

Order	Extraction	L:G Ratio	Mass of	% of Total	R _c ¹	Weight average	Number average	PDI ²
	Solvent	of PLGA ¹	PLGA (mg)	PLGA		Mol Wt (Da) ²	Mol Wt $(Da)^2$	
1	Xylenes	84.0:16.0	9.3	6.3	0.342	$17,552 \pm 333$	$15,616 \pm 266$	1.124 ± 0.002
2	Isopentyl acetate	82.9:17.1	24.8	16.4	0.425	$28,339 \pm 60$	$24,568 \pm 113$	1.148 ± 0.003
3	Toluene	82.9:17.1	37.3	24.6	0.464	$51,260 \pm 81$	$45,027 \pm 189$	1.138 ± 0.004
4	n-Butyl acetate	74.3:25.7	19.3	12.7	0.802	$26,121 \pm 184$	$22,690 \pm 179$	1.151 ± 0.001
5	2-Pentanone	72.6:27.4	22.0	14.5	0.874	$37,178 \pm 90$	$32,489 \pm 171$	1.144 ± 0.004
6	Butanone	70.9:29.1	38.2	25.2	1.00	$55,256 \pm 430$	$49,312 \pm 408$	1.120 ± 0.003
	Butanone residual	70.9:29.1	0.5	0.3	¹ NMR			
Total: 151.4 mg 100.0			100.0%	² GPC-4D				

- Trelstar 6-month comprised of multiple PLGA-types.
- Semi-solvent extraction enables separating out each specific type (L:G ratio) of PLGA and analyzing the specific fraction independently.
 Skidmore, Sarah, Justin Hadar, John Garner, Haesun Park, Kinam Park, Yan Wang, and Xiaohui Jiang.

Complex sameness: Separation of mixed poly (lactide-co-glycolide) s based on the lactide: glycolide ratio. Journal of Controlled Release 30 (2019) 174-184.







The work presented here supported by the FDA-OGD:

FDA. HHSF223201710123C (9/30/17 - 9/29/19)

Development of Analysis Technique for Structural Characterization of Star-Shaped Polyesters used for Drug Delivery

FDA. HHSF223201610091C (9/1/16 - 8/31/19)

Advanced Analytical Techniques for Mixed Polymer Drug-Delivery Systems

FDA. 1U01FD005168 (09/18/2014 - 08/31/2017)

Development of hydrogel-based in vitro dissolution apparatus for microparticle formulations

The content of this presentation is of the authors and does not necessarily represent the official views of the Food and Drug Administration.



www.Akinainc.com

John Garner jg@akinainc.com

3495 Kent Avenue West Lafayette, IN 47906 U.S.A. Tel: 1-765-464-0501