

White-Paper: Thermogels Additive effects and dissolution speed

August 28, 2017

Introduction

Polymers which have the correct architecture including poly(ethylene glycol) and hydrophobic polyester blocks at the right molecular weights have the capability to dissolve in low-temperature water and then turn into gels as the water is heated. This thermogelation property is related to the selective dehydration of the PEG chain at higher temperatures which cause the polymer micelles to bridge together forming the gel. This process involves hydrogen bonding and water structure within the solution and, as such, it can be affected by additives.

A functionality issue for thermogels is the excessive time required for their dissolution (typically days). This is due to the fact that, by definition, they are on the edge of being water soluble. Naturally, increasing the PEG block size can improve their water solubility but then such polymers are no longer thermogels but simply water soluble within the meaningful temperature range.

Purpose

In this research, we explore strategies to improve thermogel dissolution time as well as explore the effects of various additives on dissolution speed and gelation. We focus our efforts on known thermogel PLCL-PEG-PLCL at a constant solution concentration of 20% w/v water.

Materials and Methods

PLCL-PEG-PLCL (~1700-1500-1700 Da, 60:40 CL:LA) PolyVivo AK109 (Lot#60301ELH-A)

PEG 400Da (Fluka)

Dichloromethane (Malinckrodt)

Mannitol (Mannogen EZ Spray dried)

N-methylpyrrolidinone (NMP, Aldrich)

Dimethylsulfoxide (DMSO, Aldrich)

mPEG-P(DL)La (2000-2200 Da) PolyVivo AK009 (Lot# 70207STR-C)

Methods

Aliquots of known thermogel PLCL-PEG-PLCL (AK109, (200-300 mg)) were put into 5 ml vials containing 12 mm magnetic stir-bars along with the additives indicated in **Table 1**. These were then prepared as shown in **Table 1**.

Table 1. PLCL-PEG-PLCL aliquot preparations

Formulation	PLCL-PEG-PLCL (mg)	Additive (mg)	Treatment/Preparation
SGT170726A	238.2	None (control)	None
SGT170726B	300.7	Mannitol (309.9 mg)	Stirred under argon at 40 °C to mix and chopped up using metal spatula
SGT170726C	213.1	NMP (213 µL)	Stirred to mix ambient
SGT170726D	209.4	DMSO (209 µL)	Freeze dried in HarvestRight freeze-dryer

SGT170726E	282.2	DCM (282 μ L)	Dried under deep vacuum in sub-freezing conditions (dry-ice/acetone bath followed by drying in -60 °C freezer)
SGT170721F	238.6	PEG 400 (239 μ L)	Stirred to mix ambient
SGT170726A	204.7	mPEG-PLA (AK009) 60.6 mg	Stirred together at 45 °C under argon to mix
SGT170726B	218.1	mPEG-PLA (AK009) 109.7 mg	Stirred together at 45 °C under argon to mix
SGT170726C	229.3	mPEG-PLA (AK009) 215.2 mg	Stirred together at 45 °C under argon to mix
SGT170726D	224.7	PEG 400Da 112 μ l	Stirred together at 45 °C under argon to mix

Additionally, a sample of 8.3mg of PLCL-PEG-PLCL was prepared in a DSC pan and analyzed using a TA instruments Q2000 modulated DSC. Under 50 ml/min argon flush, the sample was chilled to -80 °C and then heated with modulation of ± 1 °C/min with a ramp of 3 °C/min up to 200 °C.

To each of these, a volume of water equivalent to 5 times the mass was added to each vial to generate a 20% w/v solution of the PLCL-PEG-PLCL (regardless of additive mass). The solutions were then stirred/shaken at room temperature and observed for dissolution. Upon dissolution, the samples were placed in a static 37 °C incubator (Quincy Labs, Inc. Model 12-140)

Results/Discussion

Dissolution/Gelation

Table 2. Below summarizes the dissolution and gelation results. Figure 1. Shows the dry-state condition of the prototypes.

Formulation	Dry Appearance	Dissolution	Gelation
SGT170721A	Viscous gel	Solution after 2 days of agitation at room temperature	Solution at room temperature, Gelled at 37 °C
SGT170721B	White clumpy powder	Gel after 2 days of agitation at room temperature	Gelled at room temperature, solid precipitate at 37 °C
SGT170721C	Solution	Solution after 3 minutes of agitation at room temperature	No gelation up to 37 °C
SGT170721D	Viscous gel	Solution after 2 days of agitation at room temperature	Solution at room temperature, Gelled at 37 °C
SGT170721E	Viscous gel with yellow layer	Sample observed to be contaminated with oil, dissolution determined to be > 2 hours	Not tested
SGT170721F	Solution	Gel after 3 minutes of agitation at room temperature	Gelled at room temperature, solid precipitate at 37 °C
SGT170726A	Viscous gel	Solution after 1 day of agitation at room temperature	No gelation up to 37 °C
SGT170726B	Viscous gel	Solution after 1 day of agitation at room temperature	No gelation up to 37 °C
SGT170726C	Viscous gel	Dissolution between 1 - 4 days of agitation at room temperature	Gel at room temperature, liquefied at 37 °C

SGT170726D	Viscous gel	Solution after 2 hours of agitation at room temperature	Solution at room temperature, Gelled at 37 °C
------------	-------------	---	---




Figure 1. Dry Form (A) 170721* series (B) 170726* series

Crystalline form

One manner to improve dissolution speed is to convert a material into a porous form with improved surface area for water access. Earlier attempts at simply freeze-drying the PLCL-PEG-PLCL from water in a conventional freeze-drier (HarvestRight) from water failed to generate a porous foam as the gel collapsed on itself. Additionally, freeze-drying from DMSO solvent (SGT170726D) failed to generate a porous gel as the .

In this research, the polymer was dissolved in dichloromethane (prototype SGT170726E) and deep vacuum dried (Welch Duoseal) with cooling in a dry-ice/acetone bath. **Figure 2** shows images from this.




Figure 2. (A) deep vacuum drying in dry-ice/acetone bath, (B) crystalline-foam PLCL-PEG-PLCL freshly pulled from dry-ice/acetone bath

The gel, however, did not maintain this form. Storing the gel overnight in the freezer at -20 °C it collapsed back down into the typical, liquid-gel form.

To explore this further, DSC was performed. Figure 2. Shows the results from this.

Sample: AK109-60301ELH-A
Size: 8.3000 mg
Method: Conventional MDSC

DSC

File: C:\...\170724-AK109-60301ELH-A.001

Run Date: 24-Jul-2017 10:53
Instrument: DSC Q2000 V24.11 Build 124




Figure 2. DSC scan of PLCL-PEG-PLCL (AK109)

As determined by DSC, the glass transition for this material is around -47.45°C , which puts it well below the temperature of conventional freezer storage. Outright melting occurs at 12.59°C indicating that warming the solid polymer to room temperature will liquefy it into a gel.

Additives

As for additives for additives, the results indicate the following trends:

NMP

NMP was applied as it is a biocompatible solvent. This solvent added at a 1:1 ratio with the polymer did leave a material which was readily water soluble but the solvent itself then interfered with gelation. This may still be applied in situations where it is expected that the system will leach the solvent away leaving only the polymer.

Mannitol

Mannitol is often used as an anti-caking agent for parental formulations. It actually provides for the capability to render the thermogel as a porous solid with high surface area simply by soaking up the gel and providing mechanical support. Despite this, it does not improve dissolution speed and acts to lower the gel-temperature likely by encouraging hydrogen bonding.

mPEG-P(DL)La (2000-2200Da)

The polymer selected for this application is known, by itself, to be fully water soluble with no notable thermogelation. It is also naturally a porous solid. The addition of this material did not enhance the dissolution rate of the thermogel and acted to prevent gelation at 37°C .

PEG 400

Predissolving the thermogel in PEG 400Da solvent acted to improve solubilization in water. Care must be taken, however, to not use too much PEG 400Da as it also acts to lower the gel temperature at high concentrations. PEG 400Da seems most effective when an equivalent volume that is 50% the mass of the thermogel is added (e.g. 0.5 ml PEG 400/1 g of thermogel) for the PLCL-PEG-PLCL AK109 type thermogel. This is suitable to improve dissolution speed to a few hours while not drastically lowering the gel point.